Indications for Use

The Osteotec Silicone Finger is a one-piece, sterile, flexible silicone elastomer implant of which the intended purpose is the replacement of the metacarpophalangeal (MCP) or proximal interphalangeal (PIP) joints of fingers. It is designed for the relief of pain and restoration of some function to fingers disabled by rheumatoid, degenerative or post-traumatic arthritis.

Intended User

The Silicone Finger Implants are intended for use by trained surgeons.

Description

The two intramedullary stems of the implant are joined by an integral flexible hinge which acts as a joint spacer.

The implant is not fixed into the bones and becomes stabilized by the encapsulation process. The implant is available in eleven sizes to meet various operative requirements and is supplied sterile.

A set of dedicated reusable instruments is available to prepare the bone for the implant and a set of reusable implant sizers is available to determine the appropriate size of implant. The sizers are supplied non-sterile and are not suitable for implantation. See Annex A for the list of devices and catalogue number

General Indications

As with any joint reconstruction procedure, the following general indications should be considered before deciding whether to use this implant.

- Good patient condition
- Adequate bone stock to support the implant
- Adequate skin coverage
- Patient cooperation
- Good neurovascular status
- Potentially functional musculotendinous system
- Available post-operative care

Clinical Indications

This implant is intended primarily for use in resection arthroplasty of the MCP joints and, in some cases, the PIP joint, in the following cases - In rheumatoid or post-traumatic disabilities of the MCP joint, with:

- Fixed or stiff MCP joints
- X-ray evidence of joint destruction or subluxation
- Ulnar drift, not correctable by more conservative treatment
- Contracted intrinsic and extrinsic musculature and ligament system
- Associated stiff interphalangeal joints

In rheumatoid, degenerative or post-traumatic disabilities of the PIP joint with:

- Destroyed or subluxed joint
- Stiffened joints which could not be corrected by a joint tissue release

Contraindications

- Infection
- Un-cooperative patient
- Inadequate condition of bone, skin or neurovascular system
- Permanently damaged tendon system
- Potentially successful conservative treatment
- Young patients with open epiphyses

Warnings and Precautions

- This implant should not be used in highly active or high-load-bearing patients.
- Osteotec's Surgical Technique is only a recommendation. Osteotec do not recommend a particular surgical technique for use with this implant. It remains the responsibility of the medical professional to evaluate and use the appropriate surgical technique based on their personal clinical training and experience.
- Only handle the implant using blunt instruments in order to avoid contamination or damage to the surface of the implant.
- Do not try to reshape the implant as this can adversely affect its structural integrity.
- Inspect the packaging before use. If the packaging is not intact, do not use the implant.
- The implant is for single use only – do not attempt to re-sterilise the implant as this can affect its material properties and can reduce the life of the implant.
- The use of instruments other than the dedicated instrument set may lead to inadequate fixation of the silicone finger.

Potential Complications and Adverse Side Effects

As with any joint replacement procedure, a number of complications or adverse side effects can potentially occur, these generally include but are not limited to:

- Implant loosening or fracture may necessitate revision surgery
- Excessive patient activity may cause implant wear or failure
- Infection, pain, swelling and inflammation may occur at the implant site
- Patient may have an allergic, immunological reaction, autoimmune disorder or histological response to the implant material requiring removal of the implant
- Implant movement and wear may generate wear particles that may cause or exacerbate synovitis or the formation of bone cysts in the tissues surrounding the implant.

Concerning Magnetic Resonance Environments

The Silicone Finger Implant has not been evaluated for safety and compatibility in the MR environment. The Silicone Finger Implant has not been tested for heating or migration in the MR environment.

Cleaning of the Sizers and Instruments

The sizers and instruments are reusable and should be thoroughly cleaned according to the following parameters:

- Remove any surface contamination with a disposable cloth/paper wipe. On completion of the surgical procedure the instruments should be soaked or kept moist by wrapping them in a moist towel, to prevent them from drying.
- Soak the instruments in a neutral pH enzymatic solution for a minimum of 10 minutes.
- Scrub thoroughly using a soft bristled brush to remove any surface contamination.
- Rinse thoroughly with distilled/deionised water for at least one minute.
- Dry the instrument with a clean towel or filtered air
- Perform a visual inspection of the instruments and verify that they are clean and dry. If necessary re-clean until it is visibly clean.

Warning: Silicone parts SHOULD NOT be cleaned in an ultrasonic bath.

Sizers: These instruments need to be checked for any visible signs of damage; the hinge needs to be checked for surface wear or visible tears in the material. Damaged items must be discarded and replaced.

Sterilisation of the Sizers and Instruments

The sizers and instruments are reusable and are supplied non-sterile. The sizers and instruments must be sterilised before use. The following steps are recommended -

- Double wrap the component in a porous material such as paper, fabric or sterilisation container with filters (in the US, use FDA-cleared sterilization wrap).
- Sterilise using moist steam sterilisation in a high-vacuum autoclave at 132°C for 4 minutes with 20 minutes of drying time (AAMI ST79:2010 guidelines recommendation) or HTM 01--01 Part C (2016) recommendation: 134°C for 3 minutes
- Other sterilisation methods and cycles may also be suitable. However, individuals and hospitals are advised to validate whichever method they deem appropriate at their institution.
- Do not use ethylene oxide to sterilise the sizers as gas residue in the sizers can cause an adverse tissue reaction.

Reprocessing of single use devices

Devices labelled as 'single use only' may not perform as intended by the manufacturer if re-used. The Silicone Finger implants are labelled as single use only and may not be re-used. A device should never be re-sterilised after contact with body tissue or fluids; once used the implant should be discarded.

Use of these devices cause irreversible changes to the micro and macro structure; consequently, performance characteristics of the device will be sub-optimal if re-used. Re-use of a single use device may lead to:

- An increased risk of infection
- Failure of the device to perform as intended
- Material degradation
- Endotoxin reactions

Storage

All implants must be stored in a clean, dry environment and be protected from sunlight and extremes in temperature.

Safe Disposal of the Device

Devices labelled as single use must be segregated from other reusable devices and should not be returned to a decontamination facility for reprocessing. The Osteotec Silicone Finger is considered to be special waste and should be managed as healthcare (clinical) waste. Medical devices that are to be disposed of, must be recorded, handled, and stored in accordance with SHPN 13 Series 2 and SHTN 3 and in accordance with Health and Safety, Carriage of Dangerous Goods and Waste Regulations.

Internal Standard Operating Procedures should be followed for the safe disposal of the packaging.

Handling of the device

Osteotec recommends opening the packaging in accordance with Internal Aseptic Technique methodology as the device has been provided in a double blister pack which is sterilised by gamma irradiation.

Summary of Safety and Clinical Performance

The SSCP for the Osteotec silicone finger implant can be found on the EUDAMED public website. The SSCP has the basic UDI-DI linked within this document.

For the US - Caution: Federal law restricts this device to sale by or on the order of a physician

Osteotec Ltd.
Newbury, Berkshire, United Kingdom, RG19 6HW.
TEL: +44 (0)1202 487 885
www.osteotec.co.uk

Obelis s.a.
Boulevard Général Wahis 53
1030 Brussels, BELGIUM

Annex A – List of Devices

Osteotec Silicone Finger Implant

Catalogue Number	Description
OSTF-00	Osteotec Silicone Finger Implant Size 00
OSTF-0	Osteotec Silicone Finger Implant Size 0
OSTF-1	Osteotec Silicone Finger Implant Size 1
OSTF-2	Osteotec Silicone Finger Implant Size 2
OSTF-3	Osteotec Silicone Finger Implant Size 3
OSTF-4	Osteotec Silicone Finger Implant Size 4
OSTF-5	Osteotec Silicone Finger Implant Size 5
OSTF-6	Osteotec Silicone Finger Implant Size 6
OSTF-7	Osteotec Silicone Finger Implant Size 7
OSTF-8	Osteotec Silicone Finger Implant Size 8
OSTF-9	Osteotec Silicone Finger Implant Size 9

Osteotec Silicone Finger Implant Sizer

Catalogue Number	Description
OSTF-S00	Osteotec Silicone Finger Implant Sizer Size 00
OSTF-S0	Osteotec Silicone Finger Implant Sizer Size 0
OSTF-S1	Osteotec Silicone Finger Implant Sizer Size 1
OSTF-S2	Osteotec Silicone Finger Implant Sizer Size 2
OSTF-S3	Osteotec Silicone Finger Implant Sizer Size 3
OSTF-S4	Osteotec Silicone Finger Implant Sizer Size 4
OSTF-S5	Osteotec Silicone Finger Implant Sizer Size 5
OSTF-S6	Osteotec Silicone Finger Implant Sizer Size 6
OSTF-S7	Osteotec Silicone Finger Implant Sizer Size 7
OSTF-S8	Osteotec Silicone Finger Implant Sizer Size 8
OSTF-S9	Osteotec Silicone Finger Implant Sizer Size 9

Annex A – List of Devices

Osteotec Silicone Finger Instruments

Catalogue Number	Description
OSTF-INS1	Osteotec Instrument Set
OSTF-AWL	Osteotec Silicone Finger Awl
OSTF-00D	Osteotec Silicone Finger Rasp - Size 00D
OSTF-0D	Osteotec Silicone Finger Rasp - Size 0
OSTF-1D	Osteotec Silicone Finger Rasp - Size 1D
OSTF-2D	Osteotec Silicone Finger Rasp - Size 2D
OSTF-3D	Osteotec Silicone Finger Rasp - Size 3D
OSTF-4D	Osteotec Silicone Finger Rasp - Size 4D
OSTF-5D	Osteotec Silicone Finger Rasp - Size 5D
OSTF-6D	Osteotec Silicone Finger Rasp - Size 6D
OSTF-7D	Osteotec Silicone Finger Rasp - Size 7D
OSTF-8D	Osteotec Silicone Finger Rasp - Size 8D
OSTF-9D	Osteotec Silicone Finger Rasp - Size 9D
OSTF-00P	Osteotec Silicone Finger Rasp - Size 00P
OSTF-0P	Osteotec Silicone Finger Rasp - Size 0P
OSTF-1P	Osteotec Silicone Finger Rasp - Size 1P
OSTF-2P	Osteotec Silicone Finger Rasp - Size 2P
OSTF-3P	Osteotec Silicone Finger Rasp - Size 3P
OSTF-4P	Osteotec Silicone Finger Rasp - Size 4P
OSTF-5P	Osteotec Silicone Finger Rasp - Size 5P
OSTF-6P	Osteotec Silicone Finger Rasp - Size 6P
OSTF-7P	Osteotec Silicone Finger Rasp - Size 7P
OSTF-8P	Osteotec Silicone Finger Rasp - Size 8P
OSTF-9P	Osteotec Silicone Finger Rasp - Size 9P

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