Indications for Use

The Osteotec Concentric Bone Graft Instruments is a robust set of reusable instrumentation, designed to harvest a 'plug' of bone graft and insert it into another area to aid bone fusion. They are designed to aid fusion of bones to reduce pain due to arthritis and/or correct deformity or instability.

Intended User

The Osteotec Concentric Bone Graft System is intended for use by trained surgeons.

Description

The 8 tube saws, ranging from 6.5mm Outside Diameter to 17mm Outside Diameter, are sized so that the inside diameter of the larger saw will cut a dowel that will compression fit into the hole formed by the outer diameter of the next size down. The saws and the accompanying instrumentation will allow accurate removal of damaged arthritic bone and insertion of a healthy bone dowel to encourage fusion.

The kit can be used with other implants, such as spider plates, to increase the amount of bone on some patients.

To minimise soft tissue damage, two different size displacers have been included in the kit. The two parts are slotted together and then pushed through the soft tissue until the bone is reached. The central part of the displacer can then be removed, and the tube saw can be inserted. This will stop the saw teeth from 'snagging' on and damaging the soft tissue and thus helps reduce the recovery time.

A spike has been included in the kit, this is to be inserted through the body of the tube saw and is used to centralise the saw when starting the cut.

To prevent the bones from moving while they are being cut, a jig has been included. K-wires are passed through the jig and are then drilled into each bone. The tube saw can then be passed through the center and the cut made.

The graft extractors are used to remove the 'plug' of bone from the end of the tube saw after it has been removed.

The kit is supplied in a dedicated instrument tray to aid protection during transport.

The materials used are not affected by repeated sterilisations.

Product Ref	Item
OST-CBG-01	Osteotec Concentric Bone Graft - Tube Saw S - size 1
OST-CBG-02	Osteotec Concentric Bone Graft - Tube Saw S - size 2
OST-CBG-03	Osteotec Concentric Bone Graft - Tube Saw S - size 3
OST-CBG-04	Osteotec Concentric Bone Graft - Tube Saw S - size 4
OST-CBG-05	Osteotec Concentric Bone Graft - Tube Saw L - size 5
OST-CBG-06	Osteotec Concentric Bone Graft - Tube Saw L - size 6
OST-CBG-07	Osteotec Concentric Bone Graft - Tube Saw L - size 7
OST-CBG-08	Osteotec Concentric Bone Graft - Tube Saw L - size 8
OST-CBG-10	Osteotec Concentric Bone Graft - Ejector L
OST-CBG-11	Osteotec Concentric Bone Graft - Ejector S
OST-CBG-15	Osteotec Concentric Bone Graft - Ergonomic Torque Handle
OST-CBG-16	Osteotec Concentric Bone Graft - Spike S
OST-CBG-17	Osteotec Concentric Bone Graft - Spike L
OST-CBG-24	Osteotec Concentric Bone Graft - Guide S
OST-CBG-25	Osteotec Concentric Bone Graft - Guide L
OST-CBG-26	Osteotec Concentric Bone Graft - Remover Handle
OST-CBG-20	Osteotec Concentric Bone Graft - Soft Tissue Displacer Small
OST-CBG-22	Osteotec Concentric Bone Graft - Soft Tissue Displacer Large
OST-CBG-SYSTEM	CBG Complete Tray with Instruments

General Indications

The instruments in the Concentric Bone Graft System are designed to aid fusion of bones to reduce pain due to arthritis and/or correct deformity or instability. The device is indicated for use in:

- Wrist, Carpal, Metacarpal and Phalangeal fusion
- Ankle, Tarsal, Metatarsal and Phalangeal fusion
- Spine, Lumbar and Cervical fusion

The instrument can also be used in other treatments including:

- Removal of hardware fragments
- Excision bone biopsy
- Fracture management

Surgical Benefits

- The concentric saws give the dowel an accurate compressive fit to encourage rapid fusion.
- A large range of saws allows versatile use in a large range of surgical applications.
- Spike, Shaft remover handle and k-wire channels make cutting of the bone stable and precise.
- Ergonomic design and accurate instruments ensure a simple a reproducible surgical procedure.

Contraindications

No specific contraindications exist for this procedure other than those generally associated with any orthopaedic surgery:

- Physiologically or psychologically inadequate patient.
- Inadequate skin, bone and / or neurovascular status.
- Presence of infection, sepsis, osteoporosis,
- Patients who are unwilling or incapable of following postoperative care instructions are contraindicated for these devices.

Warnings and Precautions

Possible complications which may occur with any surgery include swelling at the surgical site, flap sloughing, bleeding dehiscence, hematoma, increased sensitivity and pain, bone loss, redness, and local inflammation.

- Decontamination of reusable instruments or accessories should occur immediately after completion of the surgical procedure. Excess blood or debris should be wiped off to prevent it from drying onto the surface.
- All users should be qualified personnel with documented evidence of training and competency. Training should be inclusive of current applicable guidelines and standards and hospital policies.
- Do not use metal brushes or scouring pads during manual cleaning process.
- Use cleaning agents with low foaming surfactants for manual cleaning in order to see instruments in the cleaning solution. Cleaning agents must be easily rinsed from instruments to prevent residue.
- Mineral oil or silicone lubricants should not be used on any Osteotec instruments.
- Neutral pH enzymatic and cleaning agents are recommended for cleaning reusable instruments. It is very important that alkaline cleaning agents are thoroughly neutralised and rinsed from instruments.
- Surgical instruments must be dried thoroughly to prevent

rust formation, even if manufactured from high grade stainless steel.

- All instruments must be inspected for cleanliness of surfaces, joints, proper function, and wear and tear prior to sterilisation.
- The instruments have an indefinite life span and it will be left to the user's discretion as to when the instruments have come to the end of their useful life. Typical wear should be easy to identify as the saw teeth will appear blunt. An instrument near to the end of its useful working life will be less effective than a new one but will not prevent the successful outcome of the surgery.

Potential Complications and Adverse Side Effects

A variety of complications related to surgery or the use of the Osteotec Concentric Bone Graft can occur. Some of these may be severe, affecting the outcome of the surgery. Further surgical intervention may be required to remedy these complications if they occur.

- Excessive bleeding, which may occur due to poor surgical technique
- Damage to tissue or nerves within proximity to the grafting location.
- Infection
- Pain or discomfort
- Scar formation or other problems with the surgical incision
- Side effects from anesthesia or the surgical approach
- Skin swelling or irritation

Possible adverse effects are pain, discomfort, or abnormal sensations and nerve or soft tissue damage due to surgical trauma. Nerve or soft tissue damage, necrosis of bone or bone resorption, necrosis of the tissue or inadequate healing may result from surgical trauma.

Cleaning/ Disinfection Instructions

The instruments are reusable and should be thoroughly cleaned according to the following parameters:

- Remove any surface contamination with a disposable cloth/ paper wipe. On completion of the surgical procedure the instruments should be soaked or kept moist by wrapping them in a moist towel, to prevent them from drying.
- Disassemble instruments that require disassembly for cleaning
- Soak the instruments in a neutral pH enzymatic solution for a minimum of 10 minutes.
- Scrub thoroughly using a soft bristled brush to remove any surface contamination.
- Rinse thoroughly with distilled/deionised water for at least one minute.
- Dry the instrument with a clean towel or filtered air
- Perform a visual inspection of the instruments and verify that they are clean and dry. If necessary re-clean until it is visibly clean.

Cleaning Inspection:

- Inspect all instruments prior to sterilization or storage to ensure complete removal of soil from surfaces, cannulations, movable mechanisms and lumens
- If soil is still visible, re--clean and reprocess through the washer/disinfector

Automated Cleaning

Automatic cleaning may be used however this process would need to be validated by the end user.

Inspection and Function Testing:

- Visually inspect and check for damage and wear
- Cutting edges should be free from nicks and have a continuous edge
- Jaws and teeth should align properly
- All articulated instruments should have a smooth movement without excess play
- Locking mechanisms should fasten securely and close easily
- Check long thin instruments for bending and distortion

Sterilisation

The instruments are reusable and are supplied nonsterile.

The instruments must be sterilised before use. The following steps are recommended -

- Double wrap the component in a porous material such as paper, fabric or sterilisation container with filters (in the US, use FDA-cleared sterilization wrap).
- Sterilise using moist steam sterilisation in a high-vacuum autoclave at 132°C for 4 minutes with 20 minutes of drying time (AAMI ST79:2010 guidelines recommendation) or HTM 01- 01 Part C (2016) recommendation: 134°C for 3 minutes
- Other sterilisation methods and cycles may also be suitable. However, individuals and hospitals are advised to validate whichever method they deem appropriate at their institution.

Safe Disposal of the Device

The CBG kit is a reusable medical device. When the CBG kit has reached the end of service it must be disposed of in accordance to internal medical device disposal guidelines.

Surgical Instrument Warnings

For safe effective use of any Osteotec instrument, the surgeon must be familiar with the instrument, the method of application, and the recommended surgical technique. Instrument breakage or damage, as well as tissue damage, can occur when an instrument is subjected to excessive loads, excessive speeds, dense bone, improper use or unintended use. The patient must be cautioned, preferably in writing as to the risks associated with these types of instruments. Improper use of the device can lead to soft tissue and/or bone complications.

Surgical Instrument Precautions

Previous uses may have created imperfections or blunting of the instruments. This could lead to device failure. Instruments should be inspected prior to use. Protect instruments against scratching and nicking, such stress concentrations can lead to failure.

Storage Instructions

Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product packaging for signs of tampering, or water contamination.

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Osteotec Ltd.