


MALLET FINGER

Stéphanie Gouzou, MD
University Hospital of Strasbourg - France

METAHUS



design: e. monetaux.fr / 233 Ed.10 © AREX® - 2019

CE 0459

DISTRIBUTED BY

OSTEOTEC

Manufacturer:
AREX®
3 allée du Clos Tonnerre
91120 Palaiseau - FRANCE

www.arex.fr

AREX® - 91120 PALAISEAU - FRANCE — T. +33 [0] 169 412 212 _ F. +33 [0] 169 411 110 _ info@arex.fr

AREX®
Stay the course

HK² - METAHUS

Single use External Fixator

One HK² unit:

- 2 KW of same diameter (1.2 mm or 1.5 mm or 1.8 mm)
- 1 connector to joint 2 KW

One HK² box contains 10 HK² units (non sterile):
20 KW and 10 connectors (by 2, see photo below).

• HK218 and HK218F

- HK218 : 2 KW
- HK218F : 2 KW threaded
- + 1 connector for 2 KW

Ø 1.8 mm : DISTAL RADIUS - METACARPAL BONE



• HK215 and HK215F

- HK215 : 2 KW
- HK215F : 2 KW threaded
- + 1 connector for 2 KW

Ø 1.5 mm : CARPAL, METACARPE - PHALANGE



• HK212 and HK212F

- HK212 : 2 KW
- HK212F : 2 KW threaded
- + 1 connector for 2 KW

Ø 1.2 mm : CARPAL - PHALANX



HK2 pliers [ref. HK2P] was specifically designed for the using of HK² connectors.



BENDING the KW




CRIMPING the connectors



CUTTING NEATLY


DISTAL RADIUS FRACTURE

Gérard Hoel, MD
University Hospital of St Denis - Reunion Island




5th METACARPAL BASE FRACTURE **METAHUS**

Philippe Liverneaux, Pr
University Hospital of Strasbourg - France



BENNETH'S FRACTURE **METAHUS**

Philippe Liverneaux, Pr
University Hospital of Strasbourg - France



5th METACARPAL NECK FRACTURE **METAHUS**

Philippe Liverneaux, Pr
University Hospital of Strasbourg - France

